2025/11/08 10:27 1/3 shoulder mri guide

Veterinary Shoulder MRI Guide (v1)

By: Scott Johnson MS RT (R)(MR)

Positioning


- Spine Coil
- · Side of interest down
- Wedge under sternum
- Pull side up forelimb caudally
- Pull side down forelimb cranially until humerus is roughly perpendicular to sternum
- The idea is to get the forelimb in a relatively neutral position
 - PALPATE SHOULDER JOINT
 - NOTE WHICH COILS THE JOINT IS OVER
- 34 or 345 etc etc

First steps before scanning

- 1. Register patient
- 2. Adjust patient position to be either Left or Right decubitus, depending on side of interest
- 3. Open first scout image
- 4. Change spine coil selection to match the coils you noted inside the room
- 5. APPLY ALL
- 6. When in doubt, choose more coils

Second scout


- 1. Open scout
- 2. Maximize image
- 3. Click on 'report cursor'
- 4. Scroll through your images to find the shoulder joint on a transverse or a dorsal
- 5. Place cursor within the shoulder joint
- 6. Record the coordinates
 - Something like L30 P10 S20 etc
 - Will be different every patient
- 7. Input these coordinates with letters into the corresponding start/end points
- 8. Save, prep, scan



What to look for on second scout

https://wiki.virtual-scan.com/ Printed on 2025/11/08 10:27

2025/11/08 10:27 3/3 shoulder_mri_guide



From:

https://wiki.virtual-scan.com/ - Knowledge Portal

Permanent link:

https://wiki.virtual-scan.com/doku.php?id=shoulder_mri_guide&rev=1737066647

Last update: 2025/01/16 22:30

